

Determining the  $k_a$  of an acid-base indicator

## Indicator



Red                      Blue

$$K_a = [\text{H}_3\text{O}^+] [\text{In}^-]/[\text{HIn}]$$

$$[\text{H}_3\text{O}^+] = K_a [\text{HIn}]/[\text{In}^-]$$

$$(-\log [\text{H}_3\text{O}^+]) = (-\log K_a) + (\log [\text{In}^-]/[\text{HIn}])$$

$$\text{pH} = \text{p}K_a + \log [\text{In}^-]/[\text{HIn}]$$

## Quantitative Spectroscopy

- Beer's Law

$$A_{\lambda_1} = \epsilon_{\lambda_1}bc$$

$\epsilon$  is molar absorptivity (unique for a given compound at  $\lambda_1$ )


b is path length

c concentration


## Standard Practice

- Prepare standards of known concentration
- Measure absorbance at  $\lambda_{MAX}$
- Plot A vs. concentration
- Obtain slope
- Use slope (and intercept) to determine the concentration of the analyte in the unknown

## Typical Beer's Law Plot





# UV-vis spectrometer



# Spectroscopy


- Indicators give fairly broad visible absorption spectra









# Changing pH



# Absorption measurement






## Changing pH

$$\text{pH} = \text{pKa} + \log \left[ \frac{[\text{In}^-]}{[\text{HIn}]} \right]$$

$$\text{pKa} = \log \left[ \frac{[\text{In}^-]}{[\text{HIn}]} \right] + \text{pH}$$



## Absorption measurement

Must make measurements at two wavelengths

$$A_{\lambda_1} = \varepsilon_{\lambda_1 \text{HIn}} b C_{\text{HIn}} + \varepsilon_{\lambda_1 \text{In}^-} b C_{\text{In}^-}$$

$$A_{\lambda_2} = \varepsilon_{\lambda_2 \text{HIn}} b C_{\text{HIn}} + \varepsilon_{\lambda_2 \text{In}^-} b C_{\text{In}^-}$$

$$C_T = C_{\text{HIn}} + C_{\text{In}^-}$$

## Absorption measurement

At low pH:

$$A_{\lambda 1}^{\text{low}} = \varepsilon_{\lambda 1 \text{HIn}} b C_T \quad A_{\lambda 2}^{\text{low}} = \varepsilon_{\lambda 2 \text{HIn}} b C_T$$

At high pH:

$$A_{\lambda 1}^{\text{high}} = \varepsilon_{\lambda 1 \text{In}^-} b C_T \quad A_{\lambda 2}^{\text{high}} = \varepsilon_{\lambda 2 \text{In}^-} b C_T$$

Calculate  $\varepsilon$  for HIn and In<sup>-</sup> at both  $\lambda$ .

## In the lab

- Prepare Indicator solution adjusted to the proper pH with HA/A<sup>-</sup> buffer
- Measure UV-vis spectra and identify  $\lambda_1$  and  $\lambda_2$
- Calculate  $[In^-]/[Hin]$  from UV-vis data
- Calculate pH for each buffer solution
- Report pH vs  $\log [In^-]/[Hin]$  and draw the straight line.
- Calculate pKa as intercept