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
Il solvente

Il campione

Le **cuvette**:

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I PARAMETRI

Il metodo:

The screenshot displays the UV WinLab software interface. The main window is titled "C:\UVWINLAB\METHOD\SCANLI.MSC" and contains the following sections:

- SCAN:** Start wavelength: 950.0 nm, End wavelength: 330.0 nm, Data interval: 1.0 nm, Number of cycles: 1.
- OUTPUT:** Autosave: On Off, Autoprint: On Off, Autolist: On Off List..., Ordinate max.: 1,000, Ordinate min.: 0,000, Display: Overlay.
- End of Run Application: [Empty text box]
- Method info: scan method lab inorg

At the bottom of the main window, there are three buttons: "Scan.", "Inst.", and "Sample". These three buttons are circled in red. Below the main window, there are three smaller windows: "Data Region", "ResultWi...", and "Graph1". The status bar at the bottom shows "649.2 nm -0.0427 A Slit 2,0". The Windows taskbar at the very bottom shows the start button, "Graph Server v1.60", "InstrumentServUV", "Report Builder - C...", "UV WinLab [C:\UV...", "IT", and the system clock "11.18".

I PARAMETRI

Il metodo:

The screenshot displays the UV WinLab software interface. The main window is titled "C:\UVWINLAB\METHOD\SCANLI.MSC" and contains an "INSTRUMENT" control panel. The panel includes several settings: Ordinate mode is set to "A"; Scan speed is 480 nm/min; Lamp UV is turned Off; Lamp Vis is turned On; Slit is 1.00 nm; and Lamp change is 326.0 nm. At the bottom of the interface, a status bar shows the current measurement data: 649.2 nm, -0.0424, A, and Slit 2.0. The Windows taskbar at the bottom shows the Start button and several open applications, including Graph Serv..., InstrumentS..., Report Build..., UV WinLab [...], and UWisscreen..., along with the system clock showing 11.23.

UV WinLab [C:\UVWINLAB\METHOD\SCANLI.MSC]

File View Utilities Application Data handling Window Help

Start Autozero Setup

C:\UVWINLAB\METHOD\SCANLI.MSC

INSTRUMENT

Ordinate mode : A Scan speed : 480 nm/min

Lamp UV : On Off Smooth : 0 nm

Lamp Vis : On Off

Slit : 1.00 nm

Lamp change : 326.0 nm

Scan. Inst. Sample

Data Region ResultWi... Graph1

649.2 nm -0.0424 A Slit 2.0

start Graph Serv... InstrumentS... Report Build... UV WinLab [...] UWisscreen... IT 11.23

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UV WinLab [C:\UVWINLAB\METHOD\SCANLI.MSC]

File View Utilities Application Data handling Window Help

Start Autozero Setup

Result Filename : EXP01

Calculation factor : Factor Number of samples : 6

Fill down Delete row Insert row

No.	Sample Identity	Factor	Sample Info
1	EXP01	1,0000	[Cr(H2O)4Cl2][Cl]
2	EXP02	1,0000	
3	EXP03	1,0000	
4	EXP04	1,0000	
5	EXP05	1,0000	
6	EXP06	1,0000	

Scan. Inst. Sample

Data Region ResultWi... Graph1

649.2 nm -0.0422 A Slit 2,0

start Graph Serv... InstrumentS... Report Build... UV WinLab [... UWisscreen... IT 11.25

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Cliccare su **Start**

The screenshot shows the UV WinLab software interface. The main window is titled 'C:\UVWINLAB\METHOD\SCANLI.MSC'. It features a menu bar (File, View, Utilities, Application, Data handling, Window, Help) and a toolbar with buttons for 'Start', 'Autozero', and 'Setup'. A red arrow points to the 'Start' button. Below the toolbar, there are input fields for 'Result Filename' (EXP01), 'Calculation factor' (Factor), and 'Number of samples' (6). A table with columns 'No.', 'Sample Identity', 'Factor', and 'Sample Info' is displayed, containing six rows of data. The status bar at the bottom shows '649.2 nm -0.0422 A Slit 2.0'.

No.	Sample Identity	Factor	Sample Info
1	EXP01	1,0000	[Cr(H2O)4Cl2][Cl]
2	EXP02	1,0000	
3	EXP03	1,0000	
4	EXP04	1,0000	
5	EXP05	1,0000	
6	EXP06	1,0000	

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Cliccare su **OK**

The screenshot displays the UV WinLab software interface. At the top, the title bar reads 'UV WinLab [C:\UVWINLAB\METHOD\SCANLI.MSC]'. Below it is a menu bar with 'File', 'View', 'Utilities', 'Application', 'Data handling', 'Window', and 'Help'. A toolbar contains various icons, including a red 'Stop' button, 'Autozero', and 'ABC'. The main window features a 'Graph1' plot with a grid. The y-axis ranges from 0,00 to 1,00, and the x-axis ranges from 330,0 to 700,0. A 'Sample Change' dialog box is open in the center, with the text 'Please insert next sample: 'Blank'' and two buttons: 'OK' and 'Annulla'. A red arrow points from the text 'Cliccare su OK' to the 'OK' button. Below the graph, there are tabs for 'Scan.', 'Inst.', and 'Sample'. At the bottom, a status bar shows '649.2 nm', 'A', and 'Slit 1.0'. The Windows taskbar at the very bottom shows the 'start' button and several open applications, including 'UV WinLab [...]' and 'UVWisscreen...'. The system clock shows '11.29'.

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Cliccare su **OK**

The screenshot displays the UV WinLab software interface. At the top, the title bar reads 'UV WinLab [C:\UVWINLAB\METHOD\SCANLI.MSC]'. Below it is a menu bar with 'File', 'View', 'Utilities', 'Application', 'Data handling', 'Window', and 'Help'. A toolbar contains various icons, including a red 'Stop' button, 'Autozero', and several spectral analysis icons. The main window features a 'Graph1' plot with a y-axis from 0.00 to 1.00 and an x-axis from 330.0 to 700.0. A 'Sample Change' dialog box is open, displaying the text 'Please insert next sample: 'EXP01'' and two buttons: 'OK' and 'Annulla'. A red arrow points from the text 'Cliccare su OK' to the 'OK' button. At the bottom, a status bar shows '950.0 nm' and '0.0015 A' circled in red, along with 'Slit 1,0'. The Windows taskbar at the very bottom shows the 'start' button and several open applications, including 'Graph Serv...', 'InstrumentS...', 'Report Build...', 'UV WinLab [...]', 'UVWisscreen...', and 'IT', with the system clock showing '11:31'.

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Cliccare su
Annulla

The screenshot displays the UV WinLab software interface. At the top, the title bar reads 'UV WinLab [C:\UVWINLAB\METHOD\SCANLI.MSC]'. Below it is a menu bar with 'File', 'View', 'Utilities', 'Application', 'Data handling', 'Window', and 'Help'. A toolbar contains various icons, including a red 'Stop' button and 'Autozero'. The main window features a 'Graph1' plot with a grid. The y-axis ranges from 0,00 to 1,00, and the x-axis ranges from 330,0 to 700,0. A 'Sample Change' dialog box is open in the foreground, displaying the text 'Please insert next sample: 'EXP01'' and two buttons: 'OK' and 'Annulla'. A red arrow points from the text 'Cliccare su Annulla' to the 'Annulla' button. At the bottom of the software window, there are tabs for 'Scan.', 'Inst.', and 'Sample'. Below these tabs, a status bar shows '950.0 nm', '0.0015 A', and 'Slit 1.0'. The Windows taskbar at the very bottom shows the 'start' button and several open applications, including 'Graph Serv...', 'InstrumentS...', 'Report Build...', 'UV WinLab [...]', 'UVWisscreen...', and 'IT', with the system clock showing '11.31'.

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Cliccare su **OK**

The screenshot displays the UV WinLab software interface. The main window, titled "UV WinLab [C:\UVWINLAB\METHOD\SCANLI.MSC]", features a menu bar (File, View, Utilities, Application, Data handling, Window, Help) and a toolbar with various icons. A central window titled "Graph1" shows a plot of Absorbance (A) versus Wavelength (nm). The y-axis ranges from 0,00 to 1,00, and the x-axis ranges from 330,0 to 500. A blue line represents the recorded spectrum, which is currently flat at 0,00. Below the graph, the file name "EXP01.SP" is visible. A dialog box titled "UV WinLab [C:\UVWINLAB\METHOD\SCANLI.MSC]" is overlaid on the graph, displaying the message "Method and sampler stopped" and an "OK" button. A red arrow points from the text "Cliccare su OK" to the "OK" button. At the bottom of the software window, a status bar shows "950,0 nm", "0,0012 A", and "Slit 1,0". The Windows taskbar at the bottom shows the start button, open applications (5 NTVDM.EXE, UVWisscreenshot.doc...), and the system tray with the date and time (11.34).

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

The screenshot shows the UV WinLab software interface. The main window title is "UV WinLab [C:\UVWINLAB\METHOD\SCANLI.MS...]" and the menu bar includes "File", "View", "Utilities", "Application", "Data handling", "Window", and "Help". The toolbar contains icons for "Start", "Autozero", "Setup", and "Format Graph". A red arrow points to the "Format Graph" icon. The "Format Graph" dialog box is open, showing the "Abscissa Range" (Left: 330.00, Right: 950.00) and "Ordinate Range" (Top: 1.000, Bottom: 0.000). The dialog also has "Aurorance" checkboxes, "Set Colors", "Set Grid", "OK", "Cancel", and "Help" buttons. The main window displays a graph with a blue bar at the bottom showing "950.0 nm" and "0.0017 A". The status bar at the bottom shows the Windows taskbar with the start button, taskbar icons for "5 NTVDM.EXE" and "UVWisscreenshot.doc...", and system tray icons for "IT" and "11.41".

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
The screenshot displays the UV WinLab software interface. The main window shows a graph titled 'Graph1' with the y-axis labeled 'A' (Absorbance) ranging from 0.00 to 1.00 and the x-axis ranging from 330.0 to 400.0. A 'Find peaks' dialog box is open, allowing the user to set parameters for peak detection. The dialog box includes the following fields and options:

- Dataset:** EXP01.SP
- Threshold:** 0.100
- Abscissa start:** 950.00
- Abscissa end:** 330.00
- Label:** Peak, Base, Both
- Display:** Abscissa, Ordinate, Both

Buttons for 'OK', 'Cancel', and 'Help' are also present. The status bar at the bottom shows the current wavelength as 950.0 nm, absorbance as 0.0018 A, and slit width as 1.0. The Windows taskbar at the bottom shows the start button, taskbar with '5 NTVDM.EXE' and 'UVWisscreenshot.doc...', and system tray with 'IT' and '11.43'.

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I PARAMETRI

The screenshot displays the Spectrum software interface for a measurement. The main window is titled "Spectrum (DESKTOP-MB0T4VS-Instrument1 - Operator) - [Measurement]".

Top Panel: Includes icons for Open, Save, Print, Measurement, View, Search, and Edit Printform. It also features dropdown menus for "Data Processing" (Derivative, Peak Pick) and "Easy Macro" (Add Macros...).

Menu Bar: File, Edit, View, Graph, 2Dim, Tool, Window, Help.

Left Panel (Tree View): Shows the file structure for "LastBkg.ispd", including "Active Data: LastBkg", "LastBkg", and "C:\LabSolutions\LabSolutionsR\Data\SPETTRI DSCF\provamilani". Under "provamilani", there are sub-items for Properties, Interferogram, Background Interferogram, Peak Pick Table, and Second derivative.

Control Panel: Contains buttons for "BKG Scan", "Sample Scan", and "Monitor". It also has "Load Parameters", "Save Parameters", and "Instrument" buttons. A "Help" button and a "Switch to View after scanning" checkbox are also present. Below these are checkboxes for "View background" (checked), "Monitor autoscale" (checked), and "Continuous" (unchecked). The "Continuous" setting is set to "1 times" with a "Rest time:" field.

Parameter Fields: "Filename:" is set to "C:\LabSolutions\LabSolutionsR\Data\SPETTRI DSCF\pr...". "Auto increment" is unchecked, "Option" is "QATR-10 DIAMOND". "Sample name:" and "Sample ID:" fields are empty.

Main Plot Area: A large empty grid for the spectrum. The y-axis is labeled "%T" and ranges from 95.0 to 115.0. The x-axis is labeled "cm-1" and ranges from 4000 to 400.

Right Panel (Status/Parameters): Shows "Humidity", "Lamp", and "Laser" status as "S1101SC00275". Below this are "More Data" and "Advanced Instrument" tabs. Under "Advanced Instrument", "Measurement Mode" is "%Transmittance", "Apodization" is "Happ-Genzel", "No. of Scans" is 20, "Resolution" is 4, "Range (cm-1)" is set to "Min: 400" and "Max: 4000".

Status Log (Bottom Left): A list of system events including "FTIR Measurement ready", "SetParameters OK", "Scan requested", "Measurement stopped", "Data stored as default background spectrum", and "Data Processing: Derivative, Order=2, Points=9, Origin".

Bottom Bar: Shows the Windows taskbar with the search bar "Scrivi qui per eseguire la ricerca" and the system tray displaying the time "12:15" and date "04/11/2020".

I PARAMETRI

	<input type="checkbox"/> Humidity	
	<input type="checkbox"/> Lamp	
	<input type="checkbox"/> Laser	
	<input type="checkbox"/> S1101SC00275	
	<hr/>	
	<input type="button" value="More"/>	<input type="button" value="Advanced"/>
	<input type="button" value="Data"/>	<input type="button" value="Instrument"/>
	Measurement Mode	
	<input type="button" value="%Transmittance"/>	
	Apodization	
<input type="button" value="Happ-Genzel"/>		
No. of Scans		
<input type="button" value="20"/>		
Resolution		
<input type="button" value="4"/>		
Range (cm-1)		
Min.	<input type="button" value="400"/>	
Max.	<input type="button" value="4000"/>	

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Si fa con la sola
aria.

Cliccare su **BKG**
Scan

The screenshot displays the software interface for recording a baseline. At the top, there is a menu bar with options like 'View', 'Search', and 'Edit Printform'. Below the menu bar is a toolbar with various icons. The main control area contains three green buttons: 'BKG Scan', 'Sample Scan', and 'Monitor'. A red arrow points to the 'BKG Scan' button. To the right of these buttons are controls for 'Load Parameters', 'Save Parameters', and 'Instrument'. Below these are checkboxes for 'View background' and 'Monitor autoscale', and a 'Continuous' option set to '1 times'. A 'Rest time' field is also present. The 'Filename' field is set to 'C:\LabSolutions\LabSolutionsIR\Data\SPETTRI DSCF\pr...' and the 'Option' field is set to 'QATR-10 DIAMOND'. The 'Comment' and 'Sample name' fields are empty. The 'Sample ID' field is also empty. Below the control area is a large plot area with a grid. The y-axis is labeled '%T' and ranges from 95.0 to 115.0. The x-axis ranges from 4000 to 800. The plot area is currently empty.

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