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
Il solvente

Il campione

Le **cuvette**:

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# I PARAMETRI

Il metodo:

The screenshot displays the UV WinLab software interface. The main window is titled "C:\UVWINLAB\METHOD\SCANLI.MSC" and contains several sections:

- SCAN:** Parameters for the scan range and cycles.
  - Start wavelength: 950.0 nm
  - End wavelength: 330.0 nm
  - Data interval: 1.0 nm
  - Number of cycles: 1
- OUTPUT:** Settings for data handling and display.
  - Autosave:  On  Off
  - Autoprint:  On  Off
  - Autolist:  On  Off
  - Ordinate max.: 1,000
  - Ordinate min.: 0,000
  - Display: Overlay
- End of Run Application:** An empty text field.
- Method info:** scan method lab inorg

At the bottom of the main window, there are three buttons: "Scan.", "Inst.", and "Sample". These three buttons are circled in red. Below the main window, there is a status bar showing "649.2 nm -0.0427 A Slit 2,0". The Windows taskbar at the bottom shows the Start button and several open applications: Graph Server v1.60, InstrumentServUV, Report Builder - C..., UV WinLab [C:\UV..., and IT, with the system clock showing 11.18.

# I PARAMETRI

Il metodo:

The screenshot displays the UV WinLab software interface. The main window is titled "C:\UVWINLAB\METHOD\SCANLI.MSC" and contains an "INSTRUMENT" control panel. The parameters are as follows:

- Ordinate mode : A
- Scan speed : 480 nm/min
- Lamp UV :  On  Off
- Lamp Vis :  On  Off
- Smooth : 0 nm
- Slit : 1.00 nm
- Lamp change : 326,0 nm

At the bottom of the interface, a status bar displays the current measurement data: 649.2 nm, -0.0424 A, and Slit 2.0. The Windows taskbar at the bottom shows the Start button and several open applications: Graph Serv..., InstrumentS..., Report Build..., UV WinLab [...], and UWisscreen... The system clock shows 11.23.

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UV WinLab [C:\UVWINLAB\METHOD\SCANLI.MSC]

File View Utilities Application Data handling Window Help

Start Autozero Setup

C:\UVWINLAB\METHOD\SCANLI.MSC

Result Filename : EXP01

Calculation factor : Factor

Number of samples : 6

Fill down Delete row Insert row

No.	Sample Identity	Factor	Sample Info
1	EXP01	1,0000	[Cr(H2O)4Cl2][Cl]
2	EXP02	1,0000	
3	EXP03	1,0000	
4	EXP04	1,0000	
5	EXP05	1,0000	
6	EXP06	1,0000	

Scan. Inst. Sample

Data Region ResultWi... Graph1

649.2 nm -0.0422 A Slit 2,0

start Graph Serv... InstrumentS... Report Build... UV WinLab [ ... UWisscreen... IT 11.25

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Cliccare su **Start**

The screenshot shows the UV WinLab software interface. The main window title is 'UV WinLab [C:\UVWINLAB\METHOD\SCANLI.MSC]'. The menu bar includes 'File', 'View', 'Utilities', 'Application', 'Data handling', 'Window', and 'Help'. The toolbar contains buttons for 'Start', 'Autozero', a spectrum icon, a graph icon, a wavelength icon, a flask icon, and 'Setup'. A red arrow points to the 'Start' button. Below the toolbar is a sub-window titled 'C:\UVWINLAB\METHOD\SCANLI.MSC' with the following fields: 'Result Filename : EXP01', 'Calculation factor : Factor', and 'Number of samples : 6'. Below these fields are buttons for 'Fill down', 'Delete row', and 'Insert row'. A table displays the following data:

No.	Sample Identity	Factor	Sample Info
1	EXP01	1,0000	[Cr(H2O)4Cl2][Cl]
2	EXP02	1,0000	
3	EXP03	1,0000	
4	EXP04	1,0000	
5	EXP05	1,0000	
6	EXP06	1,0000	

At the bottom of the sub-window are tabs for 'Scan.', 'Inst.', and 'Sample', with 'Sample' being the active tab. The bottom status bar shows '649.2 nm -0.0422 A Slit 2,0'. The Windows taskbar at the bottom shows the 'start' button and several open applications: 'Graph Serv...', 'InstrumentS...', 'Report Build...', 'UV WinLab [...]', 'UVWisscreen...', and 'IT'. The system clock shows '11.25'.

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
The screenshot displays the UV WinLab software interface. At the top, the title bar reads 'UV WinLab [C:\UVWINLAB\METHOD\SCANLI.MSC]'. Below it is a menu bar with 'File', 'View', 'Utilities', 'Application', 'Data handling', 'Window', and 'Help'. A toolbar contains various icons, including a red 'Stop' button, 'Autozero', and 'ABC'. The main window features a 'Graph1' plot with a grid. The y-axis ranges from 0,00 to 1,00, and the x-axis ranges from 330,0 to 700,0. A 'Sample Change' dialog box is open, displaying the text 'Please insert next sample: 'Blank'' and two buttons: 'OK' and 'Annulla'. A red arrow points from the text 'Cliccare su OK' to the 'OK' button. Below the graph, there are tabs for 'Scan.', 'Inst.', and 'Sample'. At the bottom, a blue status bar shows '649.2 nm', 'A', and 'Slit 1.0'. The Windows taskbar at the very bottom shows the 'start' button and several open applications, including 'Graph Serv...', 'InstrumentS...', 'Report Build...', 'UV WinLab [...]', 'UVWisscreen...', and 'IT', with the system clock showing '11.29'.

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
Cliccare su **OK**

The screenshot displays the UV WinLab software interface. The main window, titled "UV WinLab [C:\UVWINLAB\METHOD\SCANLI.MSC]", features a menu bar (File, View, Utilities, Application, Data handling, Window, Help) and a toolbar with various icons. A central window titled "Graph1" shows a plot of Absorbance (A) versus wavelength (nm). The y-axis ranges from 0.00 to 1.00, and the x-axis ranges from 330.0 to 500.0. A blue line represents the recorded spectrum, which is currently flat at 0.00. Below the graph, the file name "EXP01.SP" is visible. A dialog box titled "UV WinLab [C:\UVWINLAB\METHOD\SCANLI.MSC]" is overlaid on the graph, displaying the message "Method and sampler stopped" and an "OK" button. A red arrow points from the text "Cliccare su OK" to the "OK" button. At the bottom of the interface, a status bar shows "950.0 nm", "0.0012 A", and "Slit 1.0". The Windows taskbar at the bottom indicates the system is running NTVDM.EXE and UVWisscreenshot.doc...

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

The screenshot shows the UV WinLab software interface. The main window title is "UV WinLab [C:\UVWINLAB\METHOD\SCANLI.MSS]". The menu bar includes "File", "View", "Utilities", "Application", "Data handling", "Window", and "Help". The toolbar contains icons for "Start", "Autozero", "Setup", and "Format Graph". A red arrow points to the "Format Graph" icon. The "Format Graph" dialog box is open, showing the "Abscissa Range" (Left: 330.00, Right: 950.00) and "Ordinate Range" (Top: 1.000, Bottom: 0.000). The dialog also has "Aurorange" checkboxes, "Set Colors", "Set Grid", "OK", "Cancel", and "Help" buttons. The main window shows a graph with a blue bar at the bottom displaying "950.0 nm", "0.0017 A", and "Slit 1.0". The status bar at the bottom shows the Windows taskbar with the start button, taskbar icons for "5 NTVDM.EXE" and "UVWisscreenshot.doc...", and system tray icons for "IT" and "11.41".

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
The screenshot displays the UV WinLab software interface. The main window shows a graph titled 'Graph1' with the y-axis labeled 'A' (Absorbance) ranging from 0.00 to 1.00 and the x-axis showing values 330.0 and 400.0. A 'Find peaks' dialog box is open, allowing the user to set parameters for peak detection. The dialog box includes the following fields and options:

- Dataset :** EXP01.SP
- Threshold.....** [100] (The value 100 is highlighted in blue)
- Abscissa start.....** [950,00]
- Abscissa end.....** [330,00]
- Label :**  Peak  Base  Both
- Display :**  Abscissa  Ordinate  Both

Buttons for 'OK', 'Cancel', and 'Help' are visible on the right side of the dialog box. A red arrow points to the 'Find peaks' icon in the software's toolbar. The bottom status bar shows the current wavelength as 950.0 nm, absorbance as 0.0018 A, and slit width as 1.0. The Windows taskbar at the bottom shows the Start button, running processes (NTVDM.EXE, UVWisscreenshot.doc...), and the system clock (11.43).

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# I PARAMETRI

The screenshot displays the Spectrum software interface for a measurement. The main window is titled "Spectrum (DESKTOP-MB074VS-Instrument1 - Operator) - [Measurement]".

**Top Panel:** Includes icons for Open, Save, Print, Measurement, View, Search, and Edit Printform. It also features dropdown menus for "Data Processing" (Derivative, Peak Pick) and "Easy Macro" (Add Macros...).

**Menu Bar:** File, Edit, View, Graph, 2Dim, Tool, Window, Help.

**Toolbar:** Contains various icons for file operations and measurement controls.

**Left Panel (Tree View):** Shows the project structure for "LastBkg.ispd". It includes "Active Data: LastBkg", "LastBkg", and "C:\LabSolutions\LabSolutionsR\Data\SPETTRI DSCF\provamilani". Under "provamilani", there are sub-items: Properties, Interferogram, Background Interferogram, Peak Pick Table, and Second derivative.

**Control Panel:** Features buttons for "BKG Scan", "Sample Scan", and "Monitor". It also has "Load Parameters", "Save Parameters", and "Instrument" buttons. A "Help" button and a "Switch to View after scanning" checkbox are present. Below these are checkboxes for "View background" (checked), "Monitor autoscale" (checked), and "Continuous" (unchecked) with a value of "1" times. A "Rest time:" field is also visible.

**Parameter Fields:** "Filename:" is set to "C:\LabSolutions\LabSolutionsR\Data\SPETTRI DSCF\pr...". "Auto increment" is unchecked, and "Option" is set to "QATR-10 DIAMOND". "Sample name:" and "Sample ID:" fields are empty.

**Main Plot Area:** A large grid for plotting. The vertical axis is labeled "%T" and ranges from 95.0 to 115.0. The horizontal axis is labeled "cm-1" and ranges from 4000 to 400. The plot area is currently empty.

**Right Panel (Status/Parameters):** Includes a "Stop" button. A status area shows "Humidity", "Lamp", and "Laser" with the value "S1101SC00275". Below this are "More Data" and "Advanced Instrument" tabs. Under "More Data", "Measurement Mode" is set to "%Transmittance", "Apodization" is "Happ-Genzel", "No. of Scans" is 20, and "Resolution" is 4. The "Range (cm-1)" is set to "Min: 400" and "Max: 4000".

**Status Log (Bottom Left):** A scrollable list of system events. The log shows a sequence of events including "FTIR Measurement ready", "SetParameters OK", "Scan requested", "Measurement stopped", "Data stored as default background spectrum", and "Measurement LastBkg.ispd saved". The log ends with "12:13:24 04/11/2020 FTIR Measurement ready".

**Bottom Bar:** Shows the Windows taskbar with the search bar "Scrivi qui per eseguire la ricerca" and the system tray displaying the time "12:15" and date "04/11/2020".

# I PARAMETRI

	<input type="checkbox"/> Humidity	
	<input type="checkbox"/> Lamp	
	<input type="checkbox"/> Laser	
	<input type="checkbox"/> S1101SC00275	
	<hr/>	
	<input type="button" value="More"/>	<input type="button" value="Advanced"/>
	<input type="button" value="Data"/>	<input type="button" value="Instrument"/>
	Measurement Mode	
	<input type="button" value="%Transmittance"/>	
	Apodization	
<input type="button" value="Happ-Genzel"/>		
No. of Scans		
<input type="button" value="20"/>		
Resolution		
<input type="button" value="4"/>		
Range (cm-1)		
Min.	<input type="button" value="400"/>	
Max.	<input type="button" value="4000"/>	

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Si fa con la sola  
aria.

Cliccare su **BKG**  
**Scan**

The screenshot displays the software interface for recording a baseline. At the top, there is a menu bar with options like 'View', 'Search', and 'Edit Printform'. Below the menu bar is a toolbar with various icons. The main control area contains three green buttons: 'BKG Scan', 'Sample Scan', and 'Monitor'. A red arrow points to the 'BKG Scan' button. To the right of these buttons are controls for 'Load Parameters', 'Save Parameters', and 'Instrument'. Below these are checkboxes for 'View background' and 'Monitor autoscale', and a 'Continuous' checkbox with a value of '1' and a 'Rest time' field. The bottom section of the control area includes fields for 'Filename' (C:\LabSolutions\LabSolutionsIR\Data\SPETTRI DSCF\pr...), 'Comment', 'Auto increment', 'Option' (QATR-10 DIAMOND), 'Sample name', and 'Sample ID'. The main display area is a large grid with the y-axis labeled '%T' ranging from 95.0 to 115.0 and the x-axis ranging from 4000 to 800. The Windows taskbar is visible at the bottom of the screen.

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