

## Advanced Electrophysiology 12 March 2025 - 2 April 2025

1) **Download Clampfit** 11.4.1 (or previous version if you do not have Windows 10 or 11) at:

<https://support.moleculardevices.com/s/article/Axon-pCLAMP-11-Electrophysiology-Data-Acquisition-Analysis-Software-Download-Page>

google: clampfit download

(Windows only; for mac users you need parallels, available at

<https://www.parallels.com/>

14-day free trial available

### 2) **Transfer of set of data for analysis**

In the Electrophysiological data folder, you will find 10 folders named KO Cell1-5 and WT Cell1-5. They contain electrophysiological recordings from 5 WT and 5 *Itgb3* KO layer V pyramidal neurons from medial prefrontal cortex acute slice of mice. *Itgb3* is an ASD (autism spectrum disorder) gene expressing the cell-adhesion molecule  $\beta 3$  integrin localized specifically in layer V pyramidal neurons. Each folder contains 3 files that you can open with Clampfit. They are sequential recordings of excitatory synaptic currents from the same layer V pyramidal neuron:

- 1) Baseline: Baseline recordings (12 traces recorded every 20 s)
- 2) Drug: 36 traces (every 20 s) during application of MPEP (20  $\mu$ M; 2-Methyl-6-(phenylethynyl)pyridine hydrochloride), a non-competitive antagonist of the metabotropic glutamate receptor mGluR5
- 3) NBQX: 20 traces (every 20 s) during application of NBQX (10  $\mu$ M; 2,3-Dioxo-6-nitro-1,2,3,4-tetrahydrobenzo[f]quinoxaline-7-sulfonamide) a competitive AMPA and kainate receptor antagonist.

### 3) **Opening and visualization of the traces with Clampfit**

- a. Select sweeps
- b. Adjustment visualization window
- c. Brief overview of the analysis possibilities
- d. Sampling interval (20  $\mu$ s; 50 kHz)

#### **4) What do the traces show:**

- a. Voltage command: holding potential: -70 mV, passive properties, voltage-gated currents, mEPSCs
- b. Stimulus artefact, AMPAR synaptic currents (recordings done in the presence of Bicuculline (10  $\mu$ M) a competitive GABA<sub>A</sub> antagonist and D-APV (50  $\mu$ M; D-2-amino-5-phosphonovalerate), a competitive NMDA receptor antagonist.
- c. How to isolate currents in native systems (voltage command, pharmacology, KO) as opposed to heterologous expression systems.
- d. 50 Hz noise

#### **5) How and why the experiment was done**

- a. AMPAR excitatory synaptic currents (amplitude, rise time, decay time)
  - i. Does ablation of  $\beta$ 3 integrin affects AMPAR excitatory synaptic currents under basal conditions
  - ii. Does ablation of  $\beta$ 3 integrin affects the mGluR5-mediated regulation of AMPAR synaptic currents (rational: mGluR5-mediated LTD; mGluR5 theory of fragile X-syndrome)
- b. Paired pulse ratio as a proxy of synaptic vesicle release probability ( $P_r$ )
  - i. Does ablation of  $\beta$ 3 integrin affects  $P_r$  under basal conditions
  - ii. Does ablation of  $\beta$ 3 integrin affects  $P_r$  upon blockade of mGluR5
- c. Even traces: holding potential: +40 mV
  - i. Does ablation of  $\beta$ 3 integrin affects the AMPAR rectification index under basal conditions
  - ii. Does ablation of  $\beta$ 3 integrin affects the AMPAR rectification index upon blockade of mGluR5

#### **6) How to analyze the data**

- a. How to average traces
- b. How to offset traces
- c. How to calculate current amplitude
- d. How to calculate current rise time
- e. How to calculate current decay time constant: the weighted decay time constant
- f. How to calculate AMPAR PPR
- g. How to calculate AMPAR rectification index
- h. How to calculate the effect of MPEP
- i. NBQX: how and why to digitally subtract traces

**7) Install Igor Pro and GraphPad Prism 5 on your system**

**8) Using a computer routine to automate data analysis**

- Passive properties:  $R_s$   $R_{in}$  and  $C_m$

**9) Using Two-way ANOVA – Bonferroni post-test**

WT

KO

Baseline

MPEP