Laboratorio di Microbiologia

prof. Claudio G. Ametrano
claudiogennaro.ametrano@units.it

AA 2023-2024

prof. Francesca Malfatti
fmalfatti@units.it

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DPI

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
- Procedure**
1. Flood the heat-fixed smear with crystal violet for 1 min
 2. Add iodine solution for 1 min
 3. Decolorize with alcohol briefly — about 20 sec
 4. Counterstain with safranin for 1–2 min

- Result**
-
- All cells purple
- All cells remain purple
- Gram-positive cells are purple; gram-negative cells are colorless
- Gram-positive (G^+) cells are purple; gram-negative (G^-) cells are pink to red

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Protocollo:

1. Flood air-dried, heat-fixed smear of cells for **1 minute with crystal violet** staining reagent. Please note that the quality of the smear (too heavy or too light cell concentration) will affect the Gram Stain results.
2. Wash slide in a gentle and indirect stream of tap water for 2 seconds.
3. Flood slide with the mordant: **Gram's iodine**. Wait **1 minute**.
4. Wash slide in a gentle and indirect stream of tap water for 2 seconds.
5. Flood slide with **decolorizing agent**. Wait **15 seconds** or add drop by drop to slide until decolorizing agent running from the slide runs clear (see Comments and Tips section).
6. Flood slide with counterstain, **safranin**. Wait **30 seconds to 1 minute**.
7. Wash slide in a gentle and indirect stream of tap water until no color appears in the effluent and then blot dry with absorbent paper.
8. Observe the results of the staining procedure (brightfield microscope). At the completion of the Gram Stain, gram-negative bacteria will stain pink/red and gram-positive bacteria will stain blue/purple.

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